

Test Case	1904	1908	2001	Diff
1518b-2t1c-ethip4-ip4	4.0 ± 0.0	4.0 ± 0.0	4.0 ± 0.0	+0.0 ± 0.0
1518b-2t1c-eth-l2xcba	3.9 ± 0.0	3.9 ± 0.0	3.9 ± 0.0	-0.7 ± 0.4
imix-2t1c-ethip4-ip4b	16.4 ± 0.0	16.4 ± 0.0	12.6 ± 0.2	-23.0 ± 1.0
imix-2t1c-eth-l2xcbas	16.4 ± 0.0	16.4 ± 0.0	12.0 ± 0.1	-26.7 ± 0.6
64b-2t1c-ethip4-ip4ba	27.8 ± 0.0	27.8 ± 0.0	12.3 ± 0.2	-55.7 ± 0.6
64b-2t1c-eth-l2xcbase	27.8 ± 0.0	27.8 ± 0.0	11.9 ± 0.2	-57.1 ± 0.8