

Test Case	1904	1908	2001	Diff
1518b-4t2c-ethip4-ip4	4.0 ± 0.0	4.0 ± 0.0	4.0 ± 0.0	+0.0 ± 0.0
imix-4t2c-ethip4-ip4b	16.4 ± 0.0	16.4 ± 0.0	16.4 ± 0.0	-0.0 ± 0.0
imix-4t2c-eth-l2xcbas	16.3 ± 0.2	16.4 ± 0.0	16.4 ± 0.0	-0.0 ± 0.0
1518b-4t2c-eth-l2xcba	3.9 ± 0.0	3.9 ± 0.0	3.9 ± 0.0	-0.7 ± 0.3
64b-4t2c-ethip4-ip4ba	27.8 ± 0.0	27.8 ± 0.0	25.1 ± 0.4	-9.8 ± 1.4
64b-4t2c-eth-l2xcbase	27.7 ± 0.2	27.8 ± 0.0	23.3 ± 0.7	-16.3 ± 2.5