

Test Case	2001	2005	2009	Diff
imix-2t1c-eth-12xcbas	16.4 ± 0.0	16.3 ± 0.0	16.4 ± 0.0	+0.5 ± 0.2
9000b-2t1c-eth-12xcba	NT	0.7 ± 0.0	0.7 ± 0.0	+0.1 ± 0.2
1518b-2t1c-eth-12xcba	3.9 ± 0.0	3.9 ± 0.0	3.9 ± 0.0	+0.0 ± 0.0
1518b-2t1c-ethip4-ip4	3.9 ± 0.0	3.9 ± 0.0	3.9 ± 0.0	+0.0 ± 0.0
64b-2t1c-ethip4-ip4ba	37.3 ± 0.0	37.3 ± 0.0	37.3 ± 0.0	+0.0 ± 0.0
imix-2t1c-ethip4-ip4b	16.4 ± 0.0	16.4 ± 0.0	16.4 ± 0.0	+0.0 ± 0.0
9000b-2t1c-ethip4-ip4	0.7 ± 0.0	0.7 ± 0.0	0.7 ± 0.0	-0.0 ± 0.3
64b-2t1c-eth-12xcbase	37.2 ± 0.1	37.3 ± 0.0	36.9 ± 0.9	-1.2 ± 2.3