

Test Case	2001	2005	2009	Diff
imix-2t1c-ethip4-ip4b	16.4 ± 0.0	16.4 ± 0.0	16.4 ± 0.0	+0.0 ± 0.0
imix-2t1c-eth-12xcbas	16.4 ± 0.0	16.4 ± 0.0	16.4 ± 0.0	-0.1 ± 0.2
9000b-2t1c-ethip4-ip4	0.7 ± 0.0	0.7 ± 0.0	0.7 ± 0.0	-0.2 ± 0.3
9000b-2t1c-eth-12xcbas	NT	0.7 ± 0.0	0.7 ± 0.0	-0.4 ± 0.4
1518b-2t1c-ethip4-ip4	3.9 ± 0.0	3.9 ± 0.0	3.9 ± 0.0	-0.4 ± 0.5
1518b-2t1c-eth-12xcbas	3.9 ± 0.0	3.9 ± 0.0	3.9 ± 0.0	-0.5 ± 0.5
64b-2t1c-ethip4-ip4ba	37.3 ± 0.0	37.3 ± 0.0	37.1 ± 0.1	-0.6 ± 0.2
64b-2t1c-eth-12xcbases	36.7 ± 0.1	37.3 ± 0.0	28.5 ± 9.9	-23.6 ± 26.6