

Test Case	2001	2005	2009	Diff
1518b-4t2c-ethip4-ip4	4.0 ± 0.0	4.0 ± 0.0	4.0 ± 0.0	+0.0 ± 0.0
64b-4t2c-ethip4-ip4ba	25.1 ± 0.4	37.3 ± 0.0	37.3 ± 0.0	+0.0 ± 0.0
imix-4t2c-ethip4-ip4b	16.4 ± 0.0	16.4 ± 0.0	16.4 ± 0.0	+0.0 ± 0.0
1518b-4t2c-eth-12xcbba	3.9 ± 0.0	3.9 ± 0.0	3.9 ± 0.0	-0.3 ± 0.5
64b-4t2c-eth-12xcbbase	23.3 ± 0.7	14.7 ± 0.4	14.0 ± 0.2	-4.7 ± 3.0
imix-4t2c-eth-12xcbas	16.4 ± 0.0	13.8 ± 2.0	11.6 ± 0.8	-15.6 ± 13.2